



By-pass Connectors

Type A



AK-0340/1	Connector A1, reusable (light blue)
AK-0340/2	Connector A2, reusable (light green)
AK-0340/3	Connector A3, reusable (dark green)
AK-0340/4	Connector A4, reusable (white)
AK-0340/5	Connector A5, reusable (red)
AK-0340/6	Connector A6, reusable (grey)
AK-0340/7	Connector A7, reusable (dark blue)
AK-0341	Connector VIT
AK-0340/SET	Connectors set A : AK-0340/1 + AK-0340/2 + AK-0340/3 + AK-0340/4 AK-0340/5 + AK-0340/6 + AK-0340/7+ AK-0341x2 pcs


Type D



AK-0342/1	Connector D1, reusable (white)
AK-0342/2	Connector D2, reusable (green)
AK-0342/3	Connector D3, reusable (orange)
AK-0342/SET	Connectors set D : AK-0342/1 + AK-0342/2 + AK-0342/3

Tubings/Sets

Type A



AK-0384/A1	Infusion/aspiration tubing A1
AK-0384/A2	Air/fluid tubing A2
AK-0384/A3	Assisted Infusion set A3
AK-0384/A4	Tubing set A4 - AK-0384/A1 Infusion/aspiration tubing A1 - AK-0384/A2 Air/fluid tubing A2 - AK-0384/A3 Assisted infusion Set A3
AK-0384/A5	Tubing set A5 - AK-0384/A1 Infusion/aspiration tubing A1 - AK-0384/A3 Assisted infusion Set A3

Type D

AK-0384/D1	Infusion/aspiration tubing D1
AK-0384/D2	Assisted Infusion set D2
AK-0384/D3	Infusion/aspiration tubing D3
AK-0384/D4	Tubing set D4 - AK-0384/D1 Infusion/aspiration tubing D1 - AK-5291 Infusion line
AK-0384/D5	Tubing set D5 - AK-0384/D3 Infusion/aspiration tubing D3 - AK-5291 Infusion line
AK-0384/D6	Tubing set D6 - AK-0384/D2 Assisted infusion Set D2 - AK-0384/D3 Infusion/aspiration tubing D3
AK-0384/D7	Tubing set D7 - AK-0384/D2 Assisted infusion Set D2 - AK-0384/D3 Infusion/aspiration tubing D3 - AK-0384/D8 Air/fluid tubing D8
AK-0384/D8	Air/fluid tubing D8

Attention: the tubings/sets are non-specific medical devices for fluid/air connections with a maximum pressure of 20 psi. Connection to any medical devices are the sole responsibility of the user.